HYDRO

SOFT

MEDIUM

STRONG

EVRYÁL

*Ever Young with
Hyaluronic Acid*

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Hyaluronic Bio-Revitalizing and Bio-Restructuring Gel

Treatment of Skin bio-restructuring

EVRYAL Hydro is a medical device consisting of a sodium salt of hyaluronic acid in a sterile, apyrogenic, transparent viscoelastic gel. Skin aging is a natural, progressive process. It occurs in a series of physiological changes which cause loss of moisture, elasticity and luminosity of the skin. The use of the natural product, hyaluronic acid, contained in EVRYAL Hydro can restore youthfulness to the skin. EVRYAL Hydro can be used on large areas of the body for personalised anti-aging treatment aimed at prevention, restructuring and maintenance.

Complete bio-restructuring. Thanks to its composition, EVRYAL Hydro can be used to treat all types of patients on the basis of their specific needs.

SIMPLE, VERSATILE USE

Depending on a medical assessment of the patient, different injection techniques can be used: linear, lattice, micro-depot, mixed. It's indicated for the biorevitalization and bio-restructuring of the: **face, around eyes, nasogeniene, neck, décolleté and back of hands.**

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
Linear technique 3-5 sessions every 15/30 days
(depending on photoaging)

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Back of hands

Retrograde linear technique + depots 3-5 sessions every 15/30 days
(depending on photoaging)

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Cross-linked Hyaluronic Dermal Filler

A NEW GENERATION OF CROSS-LINKED HYALURONIC ACID FILLERS FOR THE TREATMENT OF SKIN BLEMISHES

EVRYAL Cross-linked is a sterile, highly tolerated, viscoelastic filler made with cross-linked hyaluronic acid, ideal for subcutaneous applications.

EVRYAL Cross-linked helps increase the volume of the dermo-epidermal tissue thanks to the natural ability of hydrophilic hyaluronic acid molecules to attract a moderate amount of water molecules and bond with them. This allows the integration of the intercellular matrix, whilst filling intradermal spaces and making the tissues firmer.

The cross-linked hyaluronic acid contained in the product makes this filling effect longer lasting and more stable. In fact, unlike linear hyaluronic acid, the skin slowly absorbs EVRYAL Cross-linked in 6 to 8 months, looking more natural for longer. In addition to safety, BDDE – the cross-linking agent – provides

effective results, thereby meeting market requirements.

EVRYAL Cross-linked is a CE 0373 injectable, apyrogenic, resorbable medical device made with cross-linked hyaluronic acid created from non-animal sourced bacterial bio-fermentation.

It is a colourless monophasic gel contained in a disposable pre-filled, ergonomic syringe with a Luer Lock adapter.

The syringe is equipped with a plastic support that makes it easy to handle. Its accurate graduation allows you to control and read the volume of injected product.

Electronic verifications carried out throughout the entire filler production and packaging process, based on the use of state of the art technology, ensure quality, purity, and safety. These production standards, together with the extremely low concentration of nucleic acids and protein, and the bacterial endotoxin rate < 0.5 EU/ml (much lower than the limit established by the European Pharmacopoeia for injectable products) minimise any hypersensitivity reaction.

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CLINICAL TRIALS

A preclinical trial conducted by ten specialists (dermatologists, plastic surgeons, and cosmetic physicians) on 145 patients (women with a mean age of 50 years) over the course of one year has shown an excellent durability of the product. The graph n. 1 shows the duration of the treatment with the three dosages. The



satisfaction questionnaires filled out by the specialists involved in the study showed an average satisfaction index (ASI) of 8/10 for both patients and physicians. The parameters taken into the consideration in the questionnaire were the malleability and tolerability of the product, its immediate effect, pain during the treatment (which resulted to be non-existent or tolerable), and the absence of oedema (graph n. 2).

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RESULTS

Three different dosages allow for the treatment of any area of the face, reducing any type of wrinkle (fine to deep), reshaping the face contour and lips, and hydrating the skin in depth.

The result is very natural and immediately visible. EVRYAL Cross-linked has an immediate lifting effect, thanks to its intrinsic viscoelastic and hydrophilic properties, which help recreate the desired volume.

EVRYAL Cross-linked is injected subcutaneously to eliminate or at least reduce both superficial and deep wrinkles and increase soft

tissue volume (lips, cheekbones, cheeks) for face countouring or reconstructive procedures (depressed scars, acne or chicken pox outcomes, facial asymmetry and atrophy, non-surgical rhinoplasty). This allows for safe, long-lasting results, considering that resorption occurs gradually over 6 to 8 months from treatment. The treatment can be safely repeated once the filling effects are significantly attenuated.

Ease of use and intrinsic safety facilitate the physician performing the treatment. After thorough cleansing and disinfection of the skin, a few microinjections in the areas to be treated is all it takes. The treatment is well tolerated and not painful.

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 Treatment of periocular wrinkles. We recommend the linear retrograde injection technique, which should trace the path of the wrinkle with extreme precision. Avoid hypercorrections. This technique requires stretching the skin between the fingers of the hand not holding the syringe, carefully avoiding surface capillaries. After the injection, it is important to massage the treated area to spread the product subcutaneously.

MEDICAL DEVICE CE 0373

CONCENTRATION	15 mg/ml
PH	6.00 - 7.50
INITIAL HYALURONIC ACID MOLECULAR WEIGHT	2,400-4,000 kDa
TOTAL ENDOTOXINS	< 0.5 EU/ml
INDICATIONS	Treatment of periocular wrinkles
INJECTION AREA	Medium-deep dermal layer, labial submucosa
CROSS-LINKING AGENT	BDDE
FREE BDDE	< 50ppb
ONE BOX CONTAINS	Sterile syringe 1x1 mL
DURATION OF RESULTS	4-6 Months
NEEDLES	2x30G 1/2"
STERILIZATION	Autoclave

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 Indicated for the treatment of perioral wrinkles and for lip and eyebrow contouring. The product is injected into the superficial dermis, avoiding hypercorrections. If necessary, the physician can perform corrections with both the serial puncture and linear retrograde techniques, according to her dexterity. We recommend massaging the treated parts.

MEDICAL DEVICE CE 0373

CONCENTRATION	20 mg/ml
PH	6.00 – 7.50
INITIAL HYALURONIC ACID MOLECULAR WEIGHT	2,400-4,000 kDa
TOTAL ENDOTOXINS	< 0.5 EU/ml
INDICATIONS	Treatment of perioral wrinkles
INJECTION AREA	Possibly the deepest dermal layer
CROSS-LINKING AGENT	BDDE
FREE BDDE	< 50ppb
ONE BOX CONTAINS	Sterile syringe 1x1 mL
DURATION OF RESULTS	6-8 Months
NEEDLES	1x27G 1/2" - 1x26G 1/2"
STERILIZATION	Autoclave

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Treatment of medium to deep face wrinkles, volume augmentation of cheekbones, chin, nose bridge, and lips. The linear retrograde technique should be chosen over the serial puncture technique. We recommend injecting the product into the deep dermal layer.

MEDICAL DEVICE CE 0373

CONCENTRATION	25 mg/ml
PH	6.00 – 7.50
INITIAL HYALURONIC ACID MOLECULAR WEIGHT	2,400-4,000 kDa
TOTAL ENDOTOXINS	< 0.5 EU/ml
INDICATIONS	Treatment of deep facial wrinkles
INJECTION AREA	Deep or supraperiosteal hypodermis
CROSS-LINKING AGENT	BDDE
FREE BDDE	< 50ppb
ONE BOX CONTAINS	Sterile syringe 1x1 mL
DURATION OF RESULTS	6-8 Months
NEEDLES	1x27G 1/2" - 1x26G 1/2"
STERILIZATION	Autoclave

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- 1) Turlier V, Delalleau A, Casas C, Rouquier A, Bianchi P, Alvarez S, Josse G, Briant A, Dahan S, Saint-Martory C, Thenius J, Bensafi-Benaouda A, Degouy A, Schmitt AM, Redoules D. Association between collagen production and mechanical stretching in dermal extracellular matrix: in vivo effect of crosslinked hyaluronic acid filler. A randomised, placebo-controlled study. *J Dermatol Sci.* 2013 Mar; 69(3):187-94.
- 2) Cohen JL, Dayan SH, Brandt FS, Nelson DB, Axford-Gately RA, Theisen MJ, Narins RS. Systematic review of clinical trials of small and large-gel-particle hyaluronic acid injectable fillers for aesthetic soft tissue augmentation. *Dermatol Surg.* 2013 Feb; 39(2):205-31.
- 3) Eccleston D, Murphy DK. Juvéderm(*) Volbella™ in the perioral area: a 12 month prospective, multicenter, open-label study. *Clin Cosmet Investig Dermatol.* 2012;5:167-72.
- 4) Ginat DT, Schatz CJ. Imaging Features of Midface Injectable Fillers and Associated Complications. *AJNR Am J Neuroradiol.* 2012 PubMed PMID: 22837310.
- 5) Kashkouli MB, Heirati A, Pakdel F, Kivash V. Diplopia after hyaluronic acid gel injection for correction of facial tear trough deformity. *Orbit.* 2012 Oct; 31(5):330-1.
- 6) Park KY, Park MK, Li K, Seo SJ, Hong CK. Combined treatment with a nonablative infrared device and hyaluronic acid filler does not have enhanced efficacy in treating nasolabial fold wrinkles. *Dermatol Surg.* 2011 Dec; 37(12):1770-5.
- 7) Carruthers A, Carruthers J, Monheit GD, Davis PG, Tardie G. Multicenter, randomized, parallel-group study of the safety and effectiveness of onabotulinumtoxin A and hyaluronic acid dermal filler (24-mg/ml smooth, cohesive gel) alone and in combination for lower facial rejuvenation. *Dermatol Surg.* 2010 Dec; 36 Suppl 4:2121-34.

BIBLIOGRAPHY USED FOR THE CLINICAL AND BIBLIOGRAPHICAL EVALUATION OF THE PRODUCT

- 1) Turlier V, Delalleau A, Casas C, Rouquier A, Bianchi P, Alvarez S, Josse G, Briant A, Dahan S, Saint-Martory C, Thenius J, Bensafi-Benaouda A, Degouy A, Schmitt AM, Redoules D. Association between collagen production and mechanical stretching in dermal extracellular matrix: in vivo effect of cross-linked hyaluronic acid filler. A randomised, placebo-controlled study. *J Dermatol Sci.* 2013 Mar; 69(3):187-94.
- 2) Cohen JL, Dayan SH, Brandt FS, Nelson DB, Axford-Gately RA, Theisen MJ, Narins RS. Systematic review of clinical trials of small and large-gel-particle hyaluronic acid injectable fillers for aesthetic soft tissue augmentation. *Dermatol Surg.* 2013 Feb; 39(2):205-31.
- 3) Eccleston D, Murphy DK. Juvéderm(*) Volbella™ in the perioral area: a 12 month prospective, multicenter, open-label study. *Clin Cosmet Investig Dermatol.* 2012;5:167-72.
- 4) Ginat DT, Schatz CJ. Imaging Features of Midface Injectable Fillers and Associated Complications. *AJNR Am J Neuroradiol.* 2012 PubMed PMID: 22837310.
- 5) Kashkouli MB, Heirati A, Pakdel F, Kivash V. Diplopia after hyaluronic acid gel injection for correction of facial tear trough deformity. *Orbit.* 2012 Oct; 31(5):330-1.
- 6) Park KY, Park MK, Li K, Seo SJ, Hong CK. Combined treatment with a nonablative infrared device and hyaluronic acid filler does not have enhanced efficacy in treating nasolabial fold wrinkles. *Dermatol Surg.* 2011 Dec; 37(12):1770-5.
- 7) Carruthers A, Carruthers J, Monheit GD, Davis PG, Tardie G. Multicenter, randomized, parallel-group study of the safety and effectiveness of onabotulinumtoxin A and hyaluronic acid dermal filler (24-mg/ml smooth, cohesive gel) alone and in combination for lower facial rejuvenation. *Dermatol Surg.* 2010 Dec; 36 Suppl 4:2121-34.

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- 8) Monheit GD, Baumann LS, Gold MH, Goldberg DJ, Goldman MP, Narins RS, Bachtell N, Garcia E, Kablik J, Gershkovich J, Burkholder D. Novel hyaluronic acid dermal filler: dermal gel extra physical properties and clinical outcomes. *Dermatol Surg.* 2010 Nov; 36 Suppl 3:1833-41.
- 9) V, Rouquier A, Black D, Josse G, Auvergnat A, Briant A, Dahan S, Gassia V, Saint-Martory C, Zakaria W, Queille-Roussel C, Grogard-Gourdon C, Thioly-Bensoussan D, Degouy A, Schmitt AM. Assessment of the clinical efficacy of a hyaluronic acid-based deep wrinkle filler using new instrumental methods. *J Cosmet Laser Ther.* 2010 Aug; 12(4):195-202.
- 10) Arlette JP, Trotter MJ. Anatomic location of hyaluronic acid filler material injected into nasolabial fold: a histologic study. *Dermatol Surg.* 2008 Jun; 34 Suppl 1:556-62; discussion 562-3.
- 11) Tezel A, Fredrickson GH. The science of hyaluronic acid dermal filler. *J Cosmet Laser Ther.* 2008 Mar; 10(1):35-42.
- 12) Bogdan Allemann I, Baumann L. Hyaluronic acid gel (Juvéderm) preparations in the treatment of facial wrinkles and folds. *Clin Interv Aging.* 2008; 3(4):629-34.
- 13) Gold MH. Use of hyaluronic acid fillers for the treatment of the aging face. *Clin Interv Aging.* 2007; 2(3):369-76.
- 14) Falcone SJ, Berg RA. Temporary polysaccharide dermal fillers: a model for persistence based on physical properties. *Dermatol Surg.* 2009 Aug; 35(8):1238-43.
- 15) Kablik J, Monheit GD, Yu L, Chang G, Gershkovich J. Comparative physical properties of hyaluronic acid dermal fillers. *Dermatol Surg.* 2009 Feb; 35 Suppl 1:302-12.

- 8) Monheit GD, Baumann LS, Gold MH, Goldberg DJ, Goldman MP, Narins RS, Bachtell N, Garcia E, Kablik J, Gershkovich J, Burkholder D. Novel hyaluronic acid dermal filler: dermal gel extra physical properties and clinical outcomes. *Dermatol Surg.* 2010 Nov; 36 Suppl 3:1833-41.
- 9) V, Rouquier A, Black D, Josse G, Auvergnat A, Briant A, Dahan S, Gassia V, Saint-Martory C, Zakaria W, Queille-Roussel C, Grogard-Gourdon C, Thioly-Bensoussan D, Degouy A, Schmitt AM. Assessment of the clinical efficacy of a hyaluronic acid-based deep wrinkle filler using new instrumental methods. *J Cosmet Laser Ther.* 2010 Aug; 12(4):195-202.
- 10) Arlette JP, Trotter MJ. Anatomic location of hyaluronic acid filler material injected into nasolabial fold: a histologic study. *Dermatol Surg.* 2008 Jun; 34 Suppl 1:556-62; discussion 562-3.
- 11) Tezel A, Fredrickson GH. The science of hyaluronic acid dermal filler. *J Cosmet Laser Ther.* 2008 Mar; 10(1):35-42.
- 12) Bogdan Allemann I, Baumann L. Hyaluronic acid gel (Juvéderm) preparations in the treatment of facial wrinkles and folds. *Clin Interv Aging.* 2008; 3(4):629-34.
- 13) Gold MH. Use of hyaluronic acid fillers for the treatment of the aging face. *Clin Interv Aging.* 2007; 2(3):369-76.
- 14) Falcone SJ, Berg RA. Temporary polysaccharide dermal fillers: a model for persistence based on physical properties. *Dermatol Surg.* 2009 Aug; 35(8):1238-43.
- 15) Kablik J, Monheit GD, Yu L, Chang G, Gershkovich J. Comparative physical properties of hyaluronic acid dermal fillers. *Dermatol Surg.* 2009 Feb; 35 Suppl 1:302-12.

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